

Avoid costly downtime with a drop-in replacement vacuum pump. Follow the steps below and save time and money!

STEP 1




Call (801) 486-1015 or visit [COREVAX.com](https://www.corevax.com) to select a qualified reconditioned vacuum pump online

A custom COREVAX container will be shipped to you with an identical reconditioned vacuum pump



STEP 2

STEP 3




Remove and complete "Declaration of Contaminants" form with information about your used vacuum pump

Cut zip tie and remove vacuum pump from plastic bag




STEP 4

STEP 5



Swap out used vacuum pump. **Note:** Your vacuum pump has been prepped with standard oil and set to 110 voltage

Return non-seized vacuum pump in container inside plastic bag and use supplied zip tie to secure



STEP 6

STAGE 7



Repack container with shipping foam and include completed "Declaration of Contaminants"

Secure lid with straps, stick pre-paid shipping label on container and contact UPS for pick up



STEP 8

STEP 9



Check "pump repair" off your maintenance list and mark your calendar to exchange your pump in the future!